

ENVOY® 500 GLUCOSE REAGENT KIT

Product no. 55345

For *in vitro* diagnostic use

CAUTION : Federal Law restricts this device to sale by or on the order of a licensed healthcare practitioner (Rx ONLY)

INTENDED USE

Envoy® 500 Glucose Reagent is for the quantitative *in vitro* diagnostic determination of glucose in human serum and plasma on Envoy 500 Series Analysers. Glucose measurements are used in the diagnosis and treatment of carbohydrate metabolism disorders including diabetes mellitus, neonatal hypoglycemia and idiopathic hypoglycemia, and of pancreatic islet cell carcinoma.

It is not intended for use in Point of Care settings.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Glucose is the main source of energy for the human body. Glucose is converted either into glycogen to be stocked in the liver or into triglycerides to be stocked in fatty tissues. Glucose concentration in blood is regulated by several hormones, including two antagonists : insulin and glucagon. Quantification of glucose in blood is used to diagnose metabolic carbohydrates disorders such as diabetes, idiopathic hypoglycemia and pancreatic disease. The main physiological troubles are linked to hyperglycaemia (type I Diabetes mellitus and type II Diabetes mellitus). Type I diabetes mellitus is insulin-dependent, and appears mainly before 30 years old. Type II diabetes mellitus is non-insulin-dependent, and usually appears after 40 years old, but can occur earlier for obese people. Other diabetes have secondary origin, and appear after endocrinal or hepatic diseases.

METHODOLOGY ^(4,5)

Enzymatic determination of glucose according to the following reactions (Trinder reaction -End point):



4-AAP = 4-Aminoantipyrine

The red quinoneimine dye absorbs at 510 nm. The final absorbance at this wavelength is proportional to the concentration of glucose in the sample.

REAGENTS

COMPOSITION

Glucose Reagent contains 13.8 mmol/L Phosphate buffer, pH 7.4; 10 mmol/L Phenol; 0.3 mmol/L 4-Aminoantipyrine; $\geq 10,000$ U/L Glucose oxidase (*Aspergillus* sp.); ≥ 700 U/L Peroxidase (horseradish); < 0.1% Sodium azide.

WARNINGS AND PRECAUTIONS

- This reagent is for professional *in vitro* diagnostic use only.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- The reagent R contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagent always flush with copious amounts of water to prevent azide buildup.
- Dispose of contents in accordance with all local, state and federal regulations.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

PREPARATION

Glucose Reagent is ready for use on the Envoy 500 Analyzer as packaged.

STORAGE AND STABILITY

Store this reagent at 2 to 8 °C and protected from light. Do not freeze. Unopened reagents are stable to the expiration dates on the bottle labels.

The reagent is stable for 28 days onboard the Envoy 500 Analyzer.

SPECIMENS ⁽¹⁾

COLLECTION AND STORAGE

- Fresh unhemolyzed serum or lithium heparinized plasma is the preferred specimen. Do not analyze whole blood.
- Venipuncture should be performed prior to the administration of drugs. Of particular note, venipuncture performed during an acetaminophen overdose situation, when, N-acetyl-p-benzoquinone imine (NAPQI) an atypical metabolic breakdown product of acetaminophen, may be present, may lead to erroneously low glucose results. Venipuncture performed during or immediately after administration of N-acetylcysteine (NAC), a drug used to treat acetaminophen overdose, or Metamizole may lead to erroneously low glucose results.
- Samples must be separated from clot or cells promptly after collection to minimize loss of glucose through glycolysis (decrease of 5-7% in one hour in whole blood at room temperature).
- For best results, use only fresh specimens. Glucose in serum is stable for 8 hours at 25 °C or 72 hours at 2 to 8 °C.

PROCEDURE

MATERIALS PROVIDED

The Envoy 500 Glucose Reagent Kit includes the following components:

8 x 49.3 mL boats of Envoy 500 Glucose Reagent.

MATERIALS REQUIRED BUT NOT PROVIDED

- Envoy 500 Serum Calibrator (product no. 55111).
- Envoy 500 Serum Controls (product no. 55131).
- Normal saline.
- Analyzer specific consumables
- General Laboratory Equipment.

REAGENT INSTALLATION AND USE

Program the instrument using the application parameters and programming instructions provided at the end of this Instructions For Use.

Refer to the Operator Manual for additional information on installing reagents and programming the analyzer, and running samples, calibrators and controls.

The Envoy 500 Glucose Reagent is ready to use as packaged.

Do not remove the caps from the bottles until you are ready to install the reagent on the analyzer. Before installing, mix the reagents by gently inverting the reagent boats several times. Record the installation date on the label and insert the unit into the designated position on the reagent tray. Let the reagent equilibrate on the instrument for at least 30 minutes before use.

CALIBRATION

Calibrate the instrument after loading new reagent lot, after maintenance and whenever quality control results fall outside established limits. Under typical use conditions, calibration factors for this test are valid for 28 days. Refer to the Operator Manual for calibration procedures.

QUALITY CONTROL

Quality control requirements should be established in accordance with local, state and/or federal regulations or accreditation requirements.

Assay at least two levels of serum control daily. Serum Controls for Envoy500 level 1 and level 2 may be used. Controls should also be assayed after performing a reagent blank, calibrating, maintaining the instrument and after loading a new reagent unit.

Controls may be assayed more frequently based on laboratory workflow and the discretion of the user.

.../...

CALCULATIONS

All calculations are performed by the instrument.
To calculate the result in SI units (mmol/L), multiply the result in conventional units (mg/dL) by 0.0555.

PERFORMANCE CHARACTERISTICS

MEASURING RANGE

Determined according to CLSI⁽⁶⁾ EP6-A protocol, the measuring range is listed below. Samples that exceed the upper limit should be diluted 1:5 with NaCl 9 g/L solution (normal saline) and re-assayed. This extended measuring range was confirmed in a study where a high concentration of glucose was spiked into native serum samples. The recovery observed did not exceed the expected recovery by $\pm 10\%$.

The «Re-run hyperactive» function performs the dilution automatically. Results take the dilution into account.

Range	Conventional Units	SI Units
Default	20.0 to 400.0 mg/dL	1.11 to 22.20 mmol/L
Hyperactive	400.0 to 2000.0 mg/dL	22.20 to 111.01 mmol/L

LIMIT OF DETECTION (LOD) AND LIMIT OF QUANTIFICATION (LOQ)

Determined according to CLSI⁽⁷⁾ EP17-A protocol, the LoD is 1.9 mg/dL (0.11 mmol/L) and LoQ is 10.0 mg/dL (0.56 mmol/L).

EXPECTED VALUES ⁽¹⁻³⁾

Published glucose reference ranges for adults are listed below. Use these ranges only as guides. Each laboratory should establish its own reference ranges.

Reference Range	Conventional Units	SI Units
Serum/plasma	74 to 106 mg/dL	4.1 to 5.9 mmol/L

PRECISION

Determined according to CLSI⁽⁸⁾ EP5-A2 protocol.

Sample	n	mean mg/dL	Within-run CV (%)	Total
Level 1	80	36.0	0.8	2.2
Level 2	80	112.8	0.5	2.3
Level 3	80	311.5	0.6	1.9

METHOD COMPARISON

A comparative study has been performed between an Envoy500 Analyzer and an FDA-approved system equipment (Glucose oxidase method) on 101 human serum samples according to CLSI⁽⁹⁾ EP9-A2 protocol. The sample concentrations were between 19.8 and 408.7 mg/dL (1.10 and 22.69 mmol/L).

The parameters of the linear regressions are as follows :

Correlation coefficient: $(r) = 1.000$

Linear regression: $y = 0.986 x - 2.0 \text{ mg/dL} (0.11 \text{ mmol/L})$

LIMITATIONS/INTERFERING SUBSTANCES

- Do not report results outside of the usable range.

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- Studies have been performed to determine the level of interference from different compounds according to CLSI⁽¹⁰⁾ EP7-A2 protocol. Recovery is within $\pm 10\%$ of initial value of glucose concentration of 36.0 mg/dL, 108.1 mg/dL and 400.0 mg/dL.

Unconjugated bilirubin: No significant interference up to 6.0 mg/dL (103 $\mu\text{mol/L}$).

Conjugated bilirubin: No significant interference up to 5.9 mg/dL (101 $\mu\text{mol/L}$).

Hemoglobin: No significant interference up to 250 mg/dL.

Triglycerides: No significant interference up to 895 mg/dL (10.11 mmol/L).

Ascorbic acid: No significant interference up to 2.0 mg/dL (114 $\mu\text{mol/L}$).

Uric acid: No significant interference up to 20.0 mg/dL (1190 $\mu\text{mol/L}$).

Methyl dopa: No significant interference up to 0.8 mg/dL (37.9 $\mu\text{mol/L}$).

L-Dopa: Induces falsely low results at therapeutic concentrations.

Tolazamide: No significant interference up to 40.0 mg/dL (1.28 mmol/L).

Acetaminophen: No significant interference up to 30.0 mg/dL (1.98 mmol/L).

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽¹¹⁾

- Results may be falsely low when the sample is taken while levels of NAC, NAPQI (a metabolite of acetaminophen (paracetamol)) or Metamizole are significant.

- Other compounds may interfere.^(12,13)

REFERENCES

- Sacks, D.B., *Carbohydrates. Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 427.
- Dods, R.F., *Diabetes Mellitus. Clinical Chemistry : Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 580 and appendix.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 268.
- Trinder, P., *Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. Ann. Clin. Biochem.*, (1969), 6, 24.
- Burrin, J.M., Price, C.P., *Measurement of blood glucose. Ann. Clin. Biochem.*, (1985), 22, 327.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI(NCCLS) document EP6-A (2003), 23 (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), 24 (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), 24 (25).
- Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), 22 (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), 25 (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature.* Acta Clin Belg., (2004), 59, 263.
- Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press., (1997).
- Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press., (1995).

Envoy is a registered trademark of ELITech Group.

GLOSSARY OF SYMBOLS

	Contents		Manufacturer		Catalog No.
	Batch Code		See instruction for use		Use by
	Date opened / Installation date		In vitro diagnostic device		Temperature Limitation
	Number of days onboard stability		Reagent		

: Modification from previous version

ENVOY® 500 GLUCOSE REAGENT KIT

APPLICATION PARAMETERS

PRIMARY PARAMETERS		CHECK PARAMETERS		SECONDARY PARAMETERS	
Code	GLU	Reagent Limit (mABS)	200	1 st Unit Serum	mg/dL
Bar-Code	Active	Curve Acceptance (%)	100	2 nd Unit Serum	Inactive
Code for Bar-Code	331	RE-RUN SERUM		1 st Unit Urine	N/A
Test Methodology	PAP	Test Limit (Conc)	400	2 nd Unit Urine	Inactive
Method	End Point	Low Test Limit (Conc)	1.9	Dynamic Blank	Inactive
Kind of Process	Linear	Initial ABS (mABS)	N/A	Needle washes	[From Settings Table]
1st Filter	510	Final ABS (mABS)	N/A	Cuvette washes	[From Settings Table]
2nd Filter	700	Max ABS Delta (mABS)	9999	Special Wash	[From Settings Table]
Reaction direction	Increasing	Prozone Check	Inactive	Instrumental Factor	1.000
REAGENTS		Normal Range	<u>Min</u> <u>Max</u>	Shift	0.000
Number of reagents	1	Man	[User defined]	Reagent Blank	Every Day
Reagent 1 Volume μL	300	Woman	[User defined]	Decimals	1
Concentrated	Inactive	Child	[User defined]	STANDARD PARAMETERS	
Reagent 2 Volume μL	N/A	Re-run hyperactive	Active	Factor	[Determined by calibration]
Concentrated	Inactive	Re-run pathological	Inactive	Minimum	200
SAMPLE		RE-RUN URINE		Maximum	500
Name	Serum Glucose	Urine	N/A	Number of Samples	1
Sample μL	3	N/A	N/A	Max Var. (%)	10
Pre-Dilution 1:	1	N/A	N/A	Timed re-run	Inactive
Post-Dilution 1:	5	N/A	N/A	N. replicates	3
TIMES		Max ABS Delta (mABS)	N/A	Reagents ABS	[Determined by Envoy]
Sample Starter	Inactive	Prozone Check	Inactive	Pos.	[From Settings Table]
Delay Time	0	Normal Range	<u>Min</u> <u>Max</u>	Conc.	[From calibrator labeling]
Reading Time	60	Man	[User defined]	ABS	[Determined by Envoy]
Reagent 1 Incubation Time	420	Woman	[User defined]	% from last calibration	100
Reagent 2 Incubation Time	N/A	Child	[User defined]		
		Re-run hyperactive	Inactive		
		Re-run pathological	Inactive		

PROGRAMMING INSTRUCTIONS

Detailed instructions for programming reagent parameters are provided in the Envoy 500 Operator Manual and Envoy500 Settings Table.

If the Envoy 500 Chemistry System is not pre-programmed, a Glucose code must first be added before the parameters can be entered. On the menu bar, select «Test → Test Directory.» A new window will open up listing all the codes for the tests that are installed on the instrument. Click on the «New Code» button, type «GLU» into the Code field and select «Save.»

To program the application parameters, check the box next to the code for the GLU test, and select the «Parameters» button located at the bottom of the window. To program standard information, click the «Standards» button located at the bottom of the window.