



October 2008

New Product Announcement Excyte Mini Automated ESR Analyzer

Vital Diagnostics, Inc. is pleased to announce the latest addition to the Excyte product line – the Excyte Mini, part no. EX-10310.

The Excyte Mini is an automated, 10 position, random access ESR analyzer, and it is perfect for low volume hospital and POL laboratories. It has the same great features as the rest of the Excyte line including 10 minute pre-indication with results in 30 minutes, random access, automatic temperature correction, and a small footprint. Also available are the optional printer and scanner.

The Excyte Mini uses the same 1.0 mL Excyte ESR tubes, in both vacuum (part no. EX-10205) and non vacuum (part no. EX-10100) varieties. This closed tube sampling and testing increases the safety of laboratory personnel by reducing exposure to blood.

The Excyte Placement Program is available for the Excyte Mini. In exchange for a three year commitment of tubes, the customer will receive the Mini at no charge. Warranty and service are included for free.

Call 1-800-345-2822 or email USsales@vitaldiagnostics.com for additional information.

27 Wellington Road
Lincoln, RI 02865

Telephone

800-345-2822
401-642-8400

Fax

401-642-9001

www.vitaldiagnostics.com